

## ■ Secura® Semi Micro Balances

신뢰할 수 있는 세미-마이크로 전자저울



- 터치스크린 디스플레이
- 수평 이탈 시, 경고 문구 및 수평 조절 가이드 표시
- USP41 약전에 따른 최적최소무게샘플량 모니터링 기능
- 신뢰할 수 없는 리스크 감지될 경우 데이터 출력 전송 차단
- isoCAL (전자동 셀프 캘리브레이션) 기능으로 항상 정확한 결과 보장
- GxP 규격을 만족하는 출력 가능
- 별도의 소프트웨어 없이 PC로 데이터 전송 가능
- 패스워드를 통한 전자저울 설정 변경 보호

Model	Readability (mg)	Weighing Capacity (g)	Stabilization time (s)	Repeatability (mg)	Linearity Deviation (mg)
<b>Secura Semi-Micro Balances</b>					
Secura26-1SKR	0.002	21	8	0.003	0.01
Secura225D-1SKR	0.01   0.01   0.1	60   120   220	6	0.02   0.04   0.07	0.1   0.1   0.2
Secura125-1SKR	0.01   0.01	60   120	6	0.02   0.04	0.1   0.1

## ■ Quintix® Semi Micro Balances

효율적이고 간편한 세미-마이크로 전자저울



- 터치스크린 디스플레이
- isoCAL(전자동 셀프 캘리브레이션) 기능으로 항상 정확한 결과 보장
- GxP 규격을 만족하는 출력 가능
- 별도의 소프트웨어 없이 PC로 데이터 전송 가능

Model	Readability (mg)	Weighing Capacity (g)	Stabilization time (s)	Repeatability (mg)	Linearity Deviation (mg)
<b>Semi-Micro Balances</b>					
Quintix125D-1SKR	0.01   0.01   0.1	40   60   120	6	0.02   0.04   0.07	0.1   0.1   0.2
Quintix65-1SKR	0.01   0.01	40   60	6	0.02   0.04	0.1   0.1
Quintix35-1SKR	0.01	30	6	0.03	0.1

### ACC for Semi Micro

YSP01SQP	Weighing pan, 80mm, slotted
YWP01SQP	Weighing pan, 90mm; includes conversion kit
YFW01SQP	Filter weighing pan, 130mm
YDK05	Density Kit for 0.01mg Balance