

VS-5600A

혐기성 미생물 작업대 | Anaerobic Workstation

제품특징

- Anaerobic Bacteriology
- 혐기 미생물의 조작 및 배양에 최적의 조건 제공
- O₂ 농도를 1% 이하로 유지
- 샘플 출입 시 자동으로 조건을 설정하는 Air Lock
- 가스 소모가 적고 유지가 간편한 압력 조절 센서
- Forced-air Incubator 내장



VS-5600A

_ Anaerobic Workstation

Optional Accessory

Petri Plate Rack	5600-05
Mixed Regulator	5600-06
N ₂ Regulator	5600-07
Mixed Gas	5600-11
N ₂ Gas	5600-12

Catalyst Pack	5600-21
Drillite Pack	5600-22
Anaerobic Indicator	5600-23
Glove	5600-92
Sleeve	5600-94



주요사양

Anaerobic Workstation VS-5600A

Dimension (in)	1000(W) × 700(D) × 750(H)mm
(out)	1445(W) × 720(D) × 1445(H)mm
Air Lock	300(W) × 300(D) × 300(H)mm
Control	Programmable Logic Control, Manual Control
Pressure Control	Automatic
Method	Positive Pressure, Hydrogen & Catalyst
Light Bank	F.L. 20W × 2EA
Standard Device	Vacuum Pump, Glove, Catalyst Pack, Drillite Pack, Vacuum Gauge, Pressure Gauge, Shelf, Incubator Petri Plate Rack 4EA
Power	220V, 5A, 60Hz

배양기(Internal)

Anaerobic Chamber 내부에 설치하는 배양기로 온도정밀도가 우수하고 Glove를 착용한 상태에서도 사용이 편리하게 제작

Internal Incubator 5600-01

Temp. Control	Digital PID Controller
Temp. Range	RT + 5°C to 70°C
Temp. Accuracy	±0.2°C at 37°C
Display	LED Display
Door	2 Door



5600-05

_ Petri Plate Rack